

AL-MAHDI

DENTAL SUPPLIES

www.al-mahdi-dental-supplies.com Email: info@al-mahdi-dental-supplies.com Tel & WhatsApp: +447770447965

PRICE LIST - 2025

*Visit our [website](#) for product details and updates. Prices are subject to change.

CE certification on all products

DENTAL CHAIRS



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/dental-chairs-1>

AY-A8000

Complete Service Option
Delivery Only Option

Free Delivery. No Import Tax. No Customs Fees.

Price:

Complete Service Option: **£7800**
(£7410 with first order discount!)

Delivery Only Option: **£4600**
(£4370 with first order discount!)

Save over 50% on Dental Chairs!

Dental Chairs - Service Options Explained

Our chairs come with two different service options: [Complete Service Option](#) or [Delivery Only Option](#).

Service options explained	Complete Service Option	Delivery Only Option
Includes shipping and delivery to your door.	INCLUDED	INCLUDED
Includes un-boxing and chair installation.	INCLUDED	X
Includes dismantling of your old chair.	INCLUDED	X
Maintenance and repairs: 1) 5 years warranty, including spare parts and free engineer services! 2) Use your local engineer for quickness. We will pay the bill with an additional 35% (subject to approval) and provide spare parts.	INCLUDED	X
Includes 5 yrs free spare parts warranty - no engineer services.	(See Above)	INCLUDED

DENTAL CHAIRS

AY-A8000 - Price: £7800

Includes delivery, installation & 5 yrs warranty! (Complete Service Option)



Product Specifications

- 1) Thermostatic water heating system
- 2) Spittoon flushing and cup filling system
- 3) Tranquil motorized folding dental chair
- 4) Tool tray with control panel (9 chair position memories), down-mounted tool holder
- 5) 90-degree rotatable ceramic spittoon
- 6) Two-volume adjustable saliva suctions (strong, weak)
- 7) Deluxe LED sensor cold light dental lamp (gradient light intensity adjustment)
- 8) Adjustable headrest
- 9) Two 3-way syringes
- 10) LED film viewer (optional)
- 11) Deluxe multifunctional pedal switch
- 12) 2 piece of water bottle
- 13) Dentist stool
- 14) Deluxe assistant control panel
- 16) Ultrasonic scaler
- 17) Upholstery: Leather
- 18) Automatic cleaning system (cleaning all dental tubes/hoses together)



INTRAORAL SCANNERS

Runyes®



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/webshop/intraoral-scanners>

LAUNCA DL-300 (Wireless)
RUNYES 3DS 3.0 (Wired)

Purchase Option
Hire Option

Free Delivery. No Import Tax. No Customs Fees.



Price:

From £9094 (£8639 with first order discount!)

NEW

Launca DL-300 Wireless

Go Wireless with Ease

- 3 batteries included
- 1 hour of continuous scanning




LAUNCA

One of the **FASTEST** and most affordable intraoral scanners on the market!

Full arch scan in less than 25s!

(See video demonstration on our website)

X-RAY PRODUCTS



LINK TO PRODUCT PAGE
<https://www.al-mahdi-dental-supplies.com/product/10797517/handy-hdr-500-hdr-600-digital-x-ray-cmos-sensor>


HANDY - HDR 500/HDR 600 - Digital X-ray (CMOS) Sensor

Free Delivery. No Import Tax. No Customs Fees.

Price:

£1289 (£1225 with first order discount!)

X-RAY PRODUCTS



LINK TO PRODUCT PAGE
<https://www.al-mahdi-dental-supplies.com/product/12049598/handy-phosphor-plate-bsp-scanner-hds-500>

HANDY - HDS 500
Phosphor Plate (PSP) Scanner



Free Delivery. No Import Tax. No Customs Fees.

Save over 70% on PSP Scanners!

Price:

£1639 (£1557 with first order discount!)

DECONTAMINATION & STERILISATION



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/webshop/decontamination-sterilisation>

Ultrasonic Cleaner - CD 4830

Free Delivery. **No** Import Tax. **No** Customs Fees.

Price:

£261

DECONTAMINATION & STERILISATION



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/10662887/high-quality-23l-class-b-getidy-dental-steriliser-autoclave>

High Quality 23L Class B Getidy Dental Steriliser/Autoclave



Free Delivery. **No** Import Tax. **No** Customs Fees.

Price:

£3896 (£3701 with first order discount!)



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/10621393/sun-18l-b-class-autoclave-inbuilt-mini-printer-usb>

Sun 18L B Class Autoclave - Inbuilt
Mini-Printer/USB

Free Delivery. **No** Import Tax. **No** Customs Fees.

Price:

£2026 (£1925 with first order discount!)

DECONTAMINATION & STERILISATION



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/11474331/4l-stainless-steel-dental-water-distiller>

4L Stainless Steel Dental Water Distiller



Price:


£209

DENTAL BURS



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/webshop/burs/diamond-burs>



Diamond Burs

Price:

£0.50 Each

Save 50% on Diamond Burs!

DENTAL BURS



*

↑ 27mm 8.8mm ↓									
	0.6 22.0	0.8 22.0	1.0 22.0	1.2 22.0	1.4 22.0	1.6 22.0	1.8 22.0	2.1 22.0	2.3 22.0
	RA 1/2 001/006	RA 1 001/008	RA 2 001/010	RA 3 001/012	RA 4 001/014	RA 5 001/016	RA 6 001/018	RA 7 001/021	RA 8 001/023

LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/10684730/slow-speed-rose-head-round-bur-tungsten-carbide-5-pack>

Slow Speed Rose Head (Round) Bur - Tungsten Carbide

Price:

£3 Each

HIGH SPEED HANDPIECES



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/12897509/x95l-red-band-1-5-increasing-speed-fiber-optic-handpiece>

X95L - Red Band 1:5 Increasing Speed Fiber Optic Handpiece

Free Delivery. No Import Tax. No Customs Fees

Price:

£262



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/11712624/vk-207qd-450-high-speed-with-quick-connector-and-led-light>

VK-207QD - 45° High Speed with Quick Connector and LED Light

Price:

£144



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/11711124/vokodak-206qd-high-speed-with-quick-connector-led-light>

VOKODAK-206QD - High Speed With Quick Connector & LED Light

Price:

£144

HIGH SPEED HANDPIECES



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/10701693/being-fb-228-fiber-optic-handpiece-with-quick-coupling>

Being - FB-228 - Fiber Optic Handpiece with Quick Coupling



Price:

£212



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/10511352/high-speed-handpiece-vk-sl-ring-led-shadowless>

VK-SL - Ring LED - Shadowless

Price:

£144

HIGH SPEED HANDPIECES



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/10666215/handpiece-set>



Handpiece Set

Price:

£164

SLOW SPEED HANDPIECES



LINK TO PRODUCT PAGE



<https://www.al-mahdi-dental-supplies.com/product/10872263/slow-speed-push-button-handpiece>

Slow Speed Push Button Handpiece

Price:

£49

SLOW SPEED HANDPIECES



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/10707244/oral-surgery-handpiece-20-degree-angle>

Oral Surgery Handpiece - 20° Angle

Price:

£93



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/10707920/dental-implant-20-1-contra-angle-handpiece-with-light>

Dental Implant 20:1 Contra-Angle
Handpiece with Light

Price:

£222

ELECTRIC MOTOR



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/13213139/being-rose-clinic-2-electric-micromotor-built-in-type>



BEING Rose Clinic 2 - Electric Micromotor - Built-in Type

Free Delivery. No Import Tax. No Customs Fees.

Price:

£491



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/product/15159081/implant-motor-c-sailor-coxo-with-20-1-contra-angle-handpiece>



Implant Motor - C-Sailor (COXO) with 20:1 Contra-Angle Handpiece

Free Delivery. No Import Tax. No Customs Fees.

Price:

£1259 (£1197 with first order discount!)

ORTHODONTICS



Direct Printed Aligners

Powered by **Graphy**



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/webshop/orthodontics>

Direct Printed Shape Memory Aligners - By Graphy: **From £597**
<https://www.al-mahdi-dental-supplies.com/graphy-direct-printed-aligners>

Thermoplastic Premium Triple Layer Aligner Sheets: **£57 (20/pack)**

Thermoplastic Sheets For Retainers: **£13.90 (20/pack)**

Thermoplastic Sheets For Whitening Trays: **£13.90 (15/pack)**

Thermoplastic Sheets For Mouthguards: **£13.90 (10/pack)**

IPR Burs: **£0.50 Each**

Aligner Foam Toothpaste: **£4.20 (£12 RRP)**
<https://www.al-mahdi-dental-supplies.com/product/12523182/aligner-foam-toothpaste>

Sold exclusively to dentists and dental clinics

See Product Page for Details

Aligner Foam toothpaste



ALIGNER FRESH

- ✓ Cleans clear aligners and retainers
- ✓ Removes aligners build up and freshens breath
- ✓ Gradually whitens teeth with no tooth sensitivity

Contains 3% hydrogen peroxide!

A more gentle whitening experience for sensitive teeth.



ENDODONTICS



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/webshop/endodontics>

Belident - HS - NiTi (Blue) Rotary Files - Pack of 5: **£14.50**

Belident Sc-PRO- NiTi Rotary Files - Pack of 6: **£14.00**

Belident - NiTi K Hand Files - Pack of 6: **£5.00**

Belident - Stainless Steel - K Hand Files - Pack of 6: **£3.90**



Gutta Percha Points - Taper 02 (for hand files) - 120pcs: **£5**

Gutta Percha Points - Taper 04 (for rotary files) - 60pcs: **£5.50**

Brushless Endo Motor - T-Fine-II (Pro): **£276**

Brushless Endo Motor - T-Fine-II (plus) - With Apex Locator: **£472**



[See Product Page for Details](#)




CONSUMABLES



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/webshop/page-2/consumables>

Microbrushes - x100



Price:

£3.95


IMPRESSION MATERIAL



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/webshop/page-2/impression-material>

Huaer - Heavy Body - Addition Type Polyvinyl Siloxen (PVS)



Price:

£29.50



Huaer - Light Body - Addition Type Polyvinyl Siloxen (PVS)

Price:

£18.50

LED CURING LIGHTS



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/webshop/led-curing-lights>

Woodpecker LED Q Curing Light

- Installed To The Dental Chair: **£102**

VOKODAK Wireless LED 1 Sec Curing Light: **£106.50**

VK-101 LED Wireless Curing Light: **£64**

[See Product Page for Details](#)



INTRAORAL CAMERAS



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/webshop/intraoral-cameras>

HDI-220C: **£587**

HDI-712D: **£860**

Free Delivery. No Import Tax. No Customs Fees.



[See Product Page for Details](#)

EXTRACTION INSTRUMENTS

Forceps - Adult Teeth



Forceps - Primary/Childrens teeth



LINK TO PRODUCT PAGE

<https://www.al-mahdi-dental-supplies.com/webshop/page-2/extraction-instruments>

Forceps - Adult Teeth
Forceps - Primary/Childrens Teeth
Forceps - 3rd Molars
Root Forceps
Elevators



3rd Molars



Root Forceps



Elevators




Price:

£14.85 Each

DENTAL INSTRUMENTS

<https://www.al-mahdi-dental-supplies.com/webshop/page-2/dental-instruments-page-1>

Implants and Measuring Instruments / Suction Tubes / Bone Gauges



Amalgam Guns & Carriers / Retainers / Bands / Strip Holders / Napkin Holders



Rubber Dam Instruments / Rubber Dam Clamps



Syringes / Tweezer Sterilising and La Instruments



Endodontic Forceps / Tissue Forceps / Haemostatic Forceps / Towel Clamps



Retractors / Mouth Gags




Orthodontic Pliers / Cutters / Rongeurs




Needle Holders / Scissors / Micro Surgery Set



Impression Trays



Hollow Wares / Instrument Trays



Dental Implants Instruments



Dental Instruments Sets




DENTINOVA - DENTAL LABORATORY

<https://www.al-mahdi-dental-supplies.com/dentinova-digital-dental-solutions>



PRICE LIST 2024



DESIGN	MILLING	PRINTING
Single Crown..... £5	Zirconium Crown..... £27	Surgical Guide £30
Crown & Abutment..... £7	Temp. PMMA Crown..... £18	Temporary Crown..... £10
Eggshell Provisional..... £7	E-Max Crown..... £32	Waxup Model £20
Fullarch Waxup..... £50	Fullarch Zirconium..... £300	Gingivectomy Guide... £25
Surgical Guide £50		Preparation Guide £25

*Price does not include shipping. Free shipping for orders over £300.

We are the UK's **First** distribution centre for **Dentinova**, one of Africa's and the Middle East's fastest growing dental laboratories!

WHY DENTINOVA?

Trusted by over 100 facilities in more than six countries for their high quality dental designs.

Dedicated team of more than 80 hand-picked, highly qualified consultants & experienced designers covering all dental subspecialties.

Providing dental designs, milling & 3D printing services 24/7/365.

Currently processing more than 15,000 designs each month while maintaining fast delivery time, less than 24 hours for designs, and 7 - 9 days for milling and 3d printing.

UNIZ - 3D PRINTER AND ACCESSORIES




**World's Fastest
Desktop Dental Powerhouse**

Product Page: <https://www.al-mahdi-dental-supplies.com/webshop/uniz-3d-printer-and-accessories>

High-Accuracy 3D Printing Materials for Digital Dentistry

Bring production in-house or improve existing workflows with industry-leading UNIZ 3D printers. Our Dental Resins empower dental labs and practices to rapidly manufacture surgical guides, splints, fixed patterns and models, clear aligner models, and full dentures.



World's 1st Direct Printed Aligner Shape Memory Aligner



World's first patented technology **Shape Memory**

Advantages over conventional (thermoformed) aligners:

- * No attachments needed in most cases (more aesthetic)
- * Variable aligner thickness
- * More comfortable for the patient
- * You can even put the aligners in boiling water!



Distributed by:

AL-MAHDI
DENTAL SUPPLIES



Holyrood
DENTAL CARE

***First* suppliers of Direct Printed Aligners (DPA)
with shape memory in the UK!**

*Certified partners of **Graphy***

MEMORYALIGNER

Powered by Graphy

MEMORYALIGNER is the UK's first clear aligner brand using Graphy's shape memory material. Visit: www.shapememoryaligners.co.uk

The only aligner that can be cleaned in boiling water.

More Hygienic



OUR TREATMENT PACKAGES

All our orthodontic treatment planning is designed by Graphy certified specialist orthodontists.

Digital Setup/Treatment Planning: £150

1. Simple Case: £590 (No refinement fees)

Up to 10 aligners per arch

2. Moderate Case: £1090 (No refinement fees)

Up to 25 aligners per arch

3. Advanced Case: £1500 (No refinement fees)

More than 25 aligners per arch

(Or less than 25 aligners but still classified as "Advanced" due to its complexity)



50% Discount On the First Order!

COURSES

<https://www.al-mahdi-dental-supplies.com/courses>

Graphy Certification Course in Edinburgh

In-House Production of Direct Printed Shape Memory Aligners



Dr Ali Kassem
Specialist Orthodontist, PhD

2 Days Hands-On Course



Edinburgh Training & Conference Venue [16 St Mary's Street, Edinburgh EH1 1SU](#)



In cooperation with Holyrood Dental Care & Al-Mahdi Dental Supplies




Booking link: <https://www.al-mahdi-dental-supplies.com/product/15246768/graphy-certification-course-in-edinburgh>

Tel: (+44)07770447965

Email: info@al-mahdi-dental-supplies.com

Online Aligner Course

LECTURER:
Dr Ali Kassem
Specialist Orthodontist, PhD
Certified by Graphy, Invisalign, Spark and Clear Correct



Suitable for:

- ◆ Dentists who want to update their knowledge in clear aligner orthodontics
- ◆ Dentists who are new to clear aligner orthodontics

Topic Outline:

- Introduction to orthodontics
- Basic biology of orthodontic tooth movement
- Orthodontic diagnosis
- Biomechanics of Aligner treatment
- Pros and cons of Aligners
- Case selection
- Practical tips

